




INSTALLATION & OPERATION MANUAL

MODEL: NEBULA NANO
CATALOG #: F-NBLA-NNO



ADA COMPLIANT

RECORD YOUR LIFT'S
SERIAL NUMBER HERE:

1

UPDATED:
1 - 2/6/2026

WWW.AQUACREEK.COM
© 2026 AQUA CREEK PRODUCTS. ALL RIGHTS RESERVED.

9889 GARRYMORE LANE
MISSOULA, MT 59808
(406) 549-0769
FAX (406) 541-8958



TABLE OF CONTENTS

SAFETY INFORMATION.....	3
ADA REQUIREMENTS FOR POOL LIFTS (US ONLY).....	4
PRODUCT OVERVIEW AND APPLICATIONS.....	5
COMPATIBLE ANCHORS	5
UNPACKING YOUR LIFT	5
PRODUCT COMPONENTS.....	6
COMPONENT REVIEW	8
PRODUCT ASSEMBLY.....	14
OPERATION INSTRUCTIONS.....	20
MAINTENANCE INSTRUCTIONS.....	22
ANNUAL INSPECTION.....	23
BASIC TROUBLESHOOTING.....	27
WARRANTY INFORMATION.....	28



IMPORTANT SAFETY INFORMATION



Safety & Installation Instructions

When installing this product, always follow basic safety precautions, including the following:

READ ALL INSTRUCTIONS BEFORE USE

WARNING – To Reduce the Risk of Injury: Children and individuals with reduced physical, sensory, or mental capabilities, or those lacking experience or knowledge, should not use this product unless closely supervised or instructed by a responsible person.

Electrical Safety:

A wire connector is provided on all bondable anchors to connect a minimum 8 AWG (8.4 mm²) solid copper conductor to an equipotential bonding system as per NEC Article 680.82.

Maintenance & Inspection:

Periodically check the appliance for signs of wear or deterioration. Ensure all nuts, bolts, and screws are tightened to maintain safe operation.

Signage Requirement:

Affix the provided 'ADA POOL/SPA LIFT' sign in a position that maintains clear and unobstructed visibility from the lift at all times, in any location where ADA compliance is required.

Battery Safety:

Remove rechargeable batteries from the appliance before charging. Install batteries with the correct polarity. If the appliance will be unused for an extended period, remove the batteries and store them on the wall-mounted charger to maintain their charge. DO NOT short-circuit the supply terminals.

REDUCE THE RISK OF DROWNING:

- Never swim alone.
- Do not use the appliance to enter an empty pool.
- Do not spray water directly onto the electronics.
- Do not allow children to play on or around the appliance. Supervise children at all times.
- DO NOT jump or dive from this appliance.

Seat may become hot; use caution when sitting. Keep the appliance covered when not in use. Inspect and clean before each use.

READ CAREFULLY BEFORE INSTALLATION:

Inspect the entire box and packing materials for all necessary parts. Before assembly, carefully read the instructions and verify all parts using the provided diagrams and parts list. Inspect all components for damage before installation. If any damage occurred during shipping, notify Aqua Creek Products LLC within three (3) days of receipt. Proper installation is critical. Improper installation will void the Aqua Creek warranty and may compromise user safety.

ADA REQUIREMENTS FOR POOL LIFTS* (US ONLY)

1009.2.1 Pool Lift Location

Pool lift shall be located where the water level does not exceed 48". If entire pool water level exceeds 48", place pool lift where convenient.

1009.2.2 Seat Location

In the raised position, the centerline of the seat shall be located over the deck a minimum of 16" from the edge of the pool.

1009.2.3 Clear Deck Space


On the side of the seat opposite the water; a clear deck space shall be provided parallel with the seat. The space shall be 36" wide minimum and shall extend forward 48" minimum from a line located 12" behind the rear edge of the seat.

1009.2.4 Seat Height

The lift shall stop at 16" - 19" measured from the deck to the top of the seat surface when in the loading position.

1009.2.8 Submerged Depth

The lift shall submerge the seat a minimum of 18" below the stationary water level



*Compliance with ADA is the responsibility of the pool owner. Visit www.ada.gov for complete guidelines.

PRODUCT OVERVIEW AND APPLICATIONS


Introducing the Nebula Nano Pool Lift, the ideal choice for slightly above-ground spas, cold plunges, and unique tubs and settings. It can clear wall heights up to 38" (98 cm) on the lowest of the adjustable chain setting. With a weight capacity of 400 lb. (181 kg), the Nebula Nano provides reliable, powered lifting for users, while its 360-degree manual rotation allows for smooth, controlled positioning.

Built with all stainless steel and a durable powder coat finish, the Nebula Nano is designed for longevity and strength. The lift is powered by a rechargeable battery for effortless lifting and comes with a wall-mount charger, making it a practical and user-friendly accessibility solution.

COMPATIBLE ANCHORING OPTIONS

The Nebula Nano will work with the following anchors. Anchor must be installed before the lift installation can begin. All anchors are sold separately. Contact your dealer or visit our website at www.aquacreek.com for additional anchor options.


STANDARD ANCHOR OPTIONS				
DESCRIPTION	ANCHOR - 1.90"X6.00" - W/KIT (FG)	ANCHOR - WOOD DECK 1.90"X6.00"	ANCHOR - 1.90"X10.00" - W/KIT	ANCHOR - SEMI-RECESSED
CATALOG NO.	F-808SA	F-842WDA	F-808SA-10	F-77SSA
PART NO.	INV01494	INV01504	INV01495	INV01486



UNPACKING YOUR LIFT

As you unpack your new lift, please inspect all parts for damage from shipping and handling as parts may have shifted in transit. Please refer to the following checklist to ensure all parts are accounted for before you begin your assembly.

PRODUCT COMPONENTS




***16** NOT SHOWN

PRODUCT OVERVIEW AND APPLICATIONS

NEBULA NANO - INVO3799					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV03647	1	1	STOCK BOX - TIMOTION	
2	INV03491	1	1	ASSEMBLY - ANCHOR SPINDLE - 1.90" DIA. SQ DRIVE	
3	INV03813	1	1	ASSEMBLY - MAINFRAME - NEBULA NANO	
4	INV03550	1	1	ASSEMBLY - SPA SERIES HARDWARE KIT	
5	INV03822	1	1	ASSEMBLY - COVER - BATTERY & CONTROL BOX - W/ SCREWS - TIMOTION (BLUE)	
6	INV03645	1	1	ASSEMBLY - CONTROL BOX MOUNT - TIMOTION	
7	INV03602	1	-	TIMOTION - TH24 - 2 BUTTON REMOTE	
8	INV03667	1	2	ACTUATOR KIT - 400MM - TIMOTION ROTATIONAL	
9	INV03801	1	1	ASSEMBLY - LIFTING ARM - 28.00"	
10	INV03802	1	1	ASSEMBLY - CONTROL ARM - 28.00"	
11	INV02883	1	1	OFFSET SWIVEL KIT - REV/SPA SERIES	
12	INV01394	1	1	WELDMENT - STANDARD SLING HANGER BAR (FG)	
13	INV01412	1	1	TURN BAR (FG)	
14	INV01383	1	1	SLING CHAIN SET (FG)	
15	INV01386	1	1	SLING W/BARS (ADULT 38"L X 20"W) (FG)	
16	INV03812	1	1	LABEL SHEET - NEBULA NANO	


COMPONENT REVIEW




STOCK BOX - TIMOTION - INV03647					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV03598	1	-	TIMOTION - TBC2 - CHARGER	
2	INV03597	1	-	BATTERY - TIMOTION 24V 5AH	
3	INV00646	1	-	LABEL #133 - PROP 65 (BATTERY)	
4	INV01804	1	-	LABEL #146 - DAILY CHARGE	
5	INV00578	1	1	SIGN - "ADA LIFT AVAILABLE" (AQUA CREEK)	
6	INV00249	1	-	DIELECTRIC GREASE PACKET - 6 GRAMS	
7	INV00303	1	-	ANTI-SEIZE - PACKET - 5 GRAMS	
8*	INV00534	1	-	CHAIR - SEAT BELT	

*CHAIR - SEAT BELT NOT APPLICABLE WITH SLING SEAT ASSEMBLY

ASSEMBLY - ANCHOR SPINDLE - 1.90" DIA. SQ DRIVE - INV03491					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV03490	1	1	WELDMENT - ANCHOR SPINDLE SQ DRIVE	
2	INV00378	1	-	ROLLER BEARING - 1.25" I.D. X 2.33" O.D. X 5/8"	
3	INV00617	1	-	NEEDLE BEARING - K50MMX55MMX20MM	
4	INV01834	1	-	RETAINING RING - EXTERNAL (SPIRAL)	




COMPONENT REVIEW




ASSEMBLY - MAINFRAME - NEBULA NANO - INV03813					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV04413	1	1	WELDMENT- MAIN FRAME - NANO SERIES (NBLA)	
2	INV00362	2	-	BC - 1/4"-20X3/4" - FLAT HEAD COUNTERSINK	
3	INV00783	1	1	REMOTE HOLDER	
4	INV03496	1	1	WELDMENT - ACTUATOR BRACKET	
5	INV00744	2	-	NN - 5/8"-11 NYLOCK NUT	
6	INV03525	2	-	BSH - 5/8"-18X6" SOCKET HEAD BOLT 18-8 SS	
7	INV00439	1	-	BRB - 3/8"X1/2"X3/16" - BRONZE BUSHING	
8	INV00461	1	-	BS - 3/8"X1 1/4"X5/16"-18 SHOULDER BOLT	
9	INV01176	1	-	WF - 5/16" FLAT WASHER	
10	INV00742	1	-	NN - 5/16"-18 NYLOCK NUT	
11	INV00980	2	-	END CAP - RECT. 2.00"X4.00"X.120"	
12	INV01835	1	-	RETAINING RING - 1-1/4" OD EXTERNAL	
13	INV03570	1	-	SPINDLE CAP FOR NEBULA (NEB)	
14	PART OF LABEL SHEET	2	1	LABEL #223 SPA SERIES - NEBULA NANO	
15	INV03495	1	2	LASER/WJ - MAST BACK PLATE	

ASSEMBLY - SPA SERIES HARDWARE KIT - INV03550					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV00450	4	-	BRB - 1/2"X5/8"X1/4" BRONZE BUSHING	
2	INV00987	2	P1	LASER/WJ - LIFTING EAR SPACER	
3	INV00986	2	1	LASER/WJ - LIFTING EAR	
4	INV00455	2	-	BS - 1/2"X2 3/4"X3/8"-16 SHOULDER BOLT	
5	INV01174	10	-	WF - 3/8" FLAT WASHER	
6	INV00739	6	-	NN - 3/8"-16 NYLOCK NUT	
7	INV00411	4	-	BH - 3/8"-16X3 1/2" HEX HEAD BOLT	




COMPONENT REVIEW


ASSEMBLY - LIFTING ARM - 28.00" - INV03801					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV03804	1	1	WELDMENT - LIFTING ARM - 28.00"	
2	INV01082	2	-	END CAP - SQR. 2.00"X.120"	
3	INV00441	4	-	BRB - 1/2"X5/8"X5/8"X1/8" BRONZE BUSHING	
4	INV00464	1	-	BS - 3/8"X2"X5/16"-18 SHOULDER BOLT	
5	INV00448	2	-	BRB - 1/2"X3/8"X3/8" BRONZE BUSHING	
6	INV01176	1	-	WF - 5/16" FLAT WASHER	
7	INV00742	1	-	NN - 5/16"-18 NYLOCK NUT	



ASSEMBLY - CONTROL ARM - 28.00" - INV03802					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV03805	1	1	WELDMENT - CONTROL ARM - 28.00"	
2	INV01082	2	-	END CAP - SQR. 2.00"X.120"	
3	INV00441	4	-	BRB - 1/2"X5/8"X5/8"X1/8" BRONZE BUSHING	




COMPONENT REVIEW




OFFSET SWIVEL KIT - REV/SPA SERIES INV02883					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV02882	1	1	WELDMENT - OFFSET HANGER - SWIVEL	
2	INV03152	1	-	NHC - 3/8"-16 HEX CASTLE NUT	
3	INV00995	1	-	RP - 1/8"X5/8" ROLL PIN	
4	INV03172	1	1	MACHINED - SWIVEL EYEBOLT	
5	INV00270	1	-	WS - 3/8"X3/8" FLAT LDPE SPACER	
6	INV00739	2	-	NN - 3/8"-16 NYLOCK NUT	
7	INV01174	3	-	WF - 3/8" FLAT WASHER	
8	INV00451	4	-	BRB - 1/2"X5/8"X3/8" BRONZE BUSHING	
9	INV00458	2	-	BS - 1/2"X3"X3/8"-16 SHOULDER BOLT	

WELDMENT - STANDARD SLING HANGER BAR (FG) - INV01394					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV01394	1	P1	WELDMENT - STANDARD SLING HANGER BAR (FG)	





COMPONENT REVIEW




SLING W/BARS (ADULT 38"L X 20"W) (FG) - INV01386					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV01389	1	P1	SLING ONLY (ADULT 38"L X 20"W) (FG)	
2	INV01385	1	P1	SLING BARS (SET OF 04) (FG)	

SLING CHAIN SET (FG) - INV01383					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV01088	4	1	CHAIN - 3/16"	
2	INV01077	8	-	QUICK LINK - 3/16"	





TURN BAR (FG) - INV01412					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV02060	1	P1	WELDMENT - TURN BAR	
2	INV01031	1	-	HAND GRIP - RUBBER 1.00" DIA.	
3	INV01170	1	-	WF - 1/4" FLAT WASHER	
4	INV01182	1	-	WL - 1/4" LOCK WASHER	
5	INV00345	1	-	BB - 1/4"-20X1/2" BUTTON HEAD BOLT	

COMPONENT REVIEW



ASSEMBLY - COVER - BATTERY & CONTROL BOX - W/ SCREWS - TIMOTION (BLUE) - INV03822					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV00471	2	-	BSH - #10-24X3/8" SOCKET HEAD BOLT	
2	INV03819	1	1	COVER - BATTERY & CONTROL BOX - TIMOTION - BLUE	

TIMOTION - TH24 - 2 BUTTON REMOTE - INV03602					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV03602	1	-	TIMOTION - TH24 - 2 BUTTON REMOTE	



ACTUATOR KIT - 400MM - TIMOTION ROTATIONAL - INV03667					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV03600	1	-	ACTUATOR - TIMOTION - 400MM	
2	INV00657	1	1	LABEL #28 - DO NOT PUSH OR PULL	
3	INV03688	1	1	LABEL #221 ACTUATOR INFO	

PRODUCT ASSEMBLY

- **READ ALL INSTRUCTIONS FULLY AND USE ALL COMPONENTS SHOWN IN EACH STEP BEFORE MOVING ON TO THE NEXT STEP**
- **ADD A SMALL AMOUNT OF ASSEMBLY LUBE TO ALL THREADS DURING ASSEMBLY (PACKET PROVIDED). THIS WILL HELP PREVENT SEIZING OF BOLTS**
- **DO NOT OVER TIGHTEN BOLTS. OVER-TIGHTENING CAN CAUSE DAMAGE AND COMPROMISE THE SAFETY OF YOUR LIFT.**
- **TOOLS REQUIRED:**


SOCKETS/WRENCHES	ALLEN KEYS/TORX	MISCELLANEOUS
1/2"	5/32"	REAMER (AS NEEDED)
9/16"	3/16"	RUBBER Mallet (AS NEEDED)
3/4"	1/4"	
	T-30	

STEP 1

INSTALL THE MAINFRAME SUBASSEMBLY INTO THE PRE-INSTALLED ANCHOR (SEE ANCHOR MANUAL FOR INSTALLATION INSTRUCTIONS). LINE UP THE HOLE ON THE SPINDLE PLATE WITH THE FIXING ANCHOR AND SECURE THE ANCHOR BOLT, INCLUDED WITH ANCHOR KIT.

STEP 2

ATTACH LIFTING BRACKETS TO THE MAST BY USING THE INCLUDED HARDWARE IN THE SPA SERIES HARDWARE KIT. BE SURE TO INCLUDE THE SPACER PLATES BETWEEN THE MOUNTING BRACKETS. INSTALL THE BRACKETS SO THAT THE "A" & "B" LABELS ARE FACING INWARD.




***USE ASSEMBLY LUBE ON ALL BOLTS**

PRODUCT ASSEMBLY

STEP 3

ATTACH THE LIFTING ARM


ATTACH THE LIFTING ARM TO THE MAINFRAME WITH THE UNUSED HARDWARE FROM STEP 2. BE SURE TO MATCH UP THE "B" STICKER ON THE ARM WITH THE "B" STICKER ON THE MAINFRAME ATTACHMENT POINT. TIGHTEN NYLOCK NUTS TO SECURE. USE THE PREINSTALLED HARDWARE TO ATTACH THE LIFTING ARM TO THE ACTUATOR.



STEP 4

ATTACH THE ACTUATOR

ATTACH THE ACTUATOR TO THE MAINFRAME USING THE HARDWARE INCLUDED IN THE MAINFRAME. ORIENTATION SHOWN BELOW.




***USE ASSEMBLY LUBE
ON ALL BOLTS**

PRODUCT ASSEMBLY

STEP 5

ATTACH THE CONTROL ARM


ATTACH THE CONTROL ARM TO THE MAINFRAME WITH THE UNUSED HARDWARE FROM STEP 2, MATCHING THE "A" STICKERS TOGETHER. TO PREVENT SCRATCHING OR CHIPPING OF THE POWDER COAT FINISH, PLACE A CLEAN, SOFT CLOTH BETWEEN THE CONTROL AND LIFTING ARMS AT THE FURTHEST POINT FROM THE MAINFRAME.



STEP 6


ATTACH OFFSET SWIVEL SEAT ATTACHMENT

USE THE INCLUDED HARDWARE TO CONNECT THE SWIVEL ATTACHMENT TO THE LIFTING AND CONTROL ARMS.



***USE ASSEMBLY LUBE
ON ALL BOLTS**

PRODUCT ASSEMBLY



STEP 7
ATTACH HANGER BAR AND SLING
SEAT ASSEMBLY

PLACE THE HANGER BAR THROUGH THE EYE BOLT. THEN

STEP 8
ATTACH SEAT TO CHAINS

ATTACH THE CHAIN LINKS TO THE MESH SEAT BY CONNECTING THE QUICK LINKS TO THE BAR INSERTS OF THE MESH SEAT (4 PLACES).

STEP 9
ATTACH CHAINS TO HANGER


LOOP THE CHAIN LINKS ONTO THE HOOKED ENDS OF THE HANGER BAR. MAKE SURE THE CHAIN LINKS ARE FREE FROM ANY KINKS.

PRODUCT ASSEMBLY

STEP 10

ATTACH TURN BAR


ATTACH TURN BAR USING THE INCLUDED HARDWARE. BE SURE TO POSITION TURN BAR TO THAT IT WILL NOT HIT THE SPA WALL, AND THE OPERATOR HAS ENOUGH CLEARANCE TO MANEUVER THE LIFT.



PRODUCT ASSEMBLY


CONGRATULATIONS!

PRODUCT ASSEMBLY IS COMPLETE! RINSE THE LIFT OFF WITH FRESH WATER AFTER ASSEMBLY IS COMPLETE TO REMOVE ANY GREASE LEFT OVER FROM THE ASSEMBLY PROCESS.




OPERATION INSTRUCTIONS


NOTE: THIS CONTROL BOX COMES WITH AN "EMERGENCY STOP" BUTTON ON TOP. LIFT WILL NOT WORK WHILE THE BUTTON IS COMPRESSED. PLEASE ENSURE BUTTON IS NOT ENGAGED WHILE LIFT IS IN OPERATION. TO DISENGAGE, PRESS AND TWIST BUTTON UNTIL IT POPS UP.



2-BUTTON REMOTE FUNCTION



ROTATIONAL UPGRADE KIT



POWERED ROTATION KIT - NEBULA - INV03576					
ITEM	PART NO.	QTY.	REV.	DESCRIPTION	REORDER
1	INV02755	1	2	ASSEMBLY - MOTOR MOUNT	
2	INV03502	1	1	ASSEMBLY - MOTOR COVER - BLUE	
3	INV03603	1	-	TIMOTION - TH24 - 4 BUTTON REMOTE	


INSTALLATION

To install the Powered Rotation Kit, remove the plastic cap on the spindle and clean out any debris that may be in the connecting square hole. Clean and lubricate the 4 mounting holes next to the spindle top. Line up the square motor shaft with the spindle and insert the motor assembly. Rotate the Nebula mainframe around until the slots on the motor bracket line up with the mounting hole on the mainframe. Insert and tighten the 4 included screws.

Remove the plastic plug in the third available port in the control box. Plug the motor's cord into this port.

Next, swap out the 2 button remote for the included 4 button remote. Ensure all connections are secure.

Protect your electronics by installing the cover. Place the cover over the motor as shown and secure with the included hardware.



MAINTENANCE INSTRUCTIONS

PROCEDURE	BEFORE EACH USE	DAILY	WEEKLY	MONTHLY
Check battery level indicator on control box.	☑			
Run the lift through a complete test cycle (fully extend actuator) to verify it is functioning properly.	☑			
Wash lift with fresh water and dry with clean, soft, non-abrasive cloth.		☑		
Recharge battery.		☑		
Visually inspect lift for damage, corrosion, and loose or missing hardware.		☑		
Check all contact points for damage and/or corrosion.			☑	
Repair, clean and apply dielectric grease to all contact points.			☑	
Thoroughly clean lift frame and apply a liberal coat of car wax to maintain the lift's finish.				☑
Check all Warning and Cautionary labels to make sure they are not faded or worn. Replace as needed.				☑

PROPER CARE NOTES:

- Use only fresh water to wash your lift. Do not wash with pool water.
- Use only clean, soft, non-abrasive cloths on your lift.
- Do not store lift in pump room or near chemicals.
- Do not allow children to play on or around this lift.
- When cleaning the lift, do not spray water directly on control box or battery.
- The lift may need to be bonded according to local code requirements.
- If the lift is not used for extended periods, remove the battery from the lift and store it on the wall charger to maintain the batteries charge. The charger should be plugged in and in a climate controlled area.

STAIN REMOVAL NOTES:

To remove stubborn stains from finish, spray affected area with mild dish soap, rub briskly with a 3M Scotch Brite™ pad. Use white epoxy paint to touch up the finish. Apply wax to protect when dry.

LIFT INSPECTION INFORMATION

Test Date:	
Pool Lift Model:	
Serial Number:	
Inspected By:	

Inspection Notes:

Please read and fill out all information carefully

Improper use and/or maintenance of pool lift can lead to actuator and/or other mechanical failures and possible personal injury. Please read and follow the product operators' manual and the instructions below.

NOTE: Completely cycle the lift without a load before each use to verify the lift is functioning properly.

- Product service life depends on the load characteristics, type and frequency of use, and the environment in which the product is being used. For example, persons not calmly sitting but falling down with their full body weight in the chair, or using the pool lift's chair as a diving platform can cause premature actuator failure and/or possible personal injury.
- Based on historical product service life analysis, Aqua Creek recommends replacing the actuator(s) and other electronic components after 5 years of service. Check component serial numbers and purchase records to verify the age of your components.
- A functional weight lifting force of a maximum of 100lbs. Must be periodically verified by applying a non-live load (without a person) through the entire operating cycle of the lift. This process is performed in order to verify the smooth operation. Make note of any unusual noises or sudden movements coming from the actuator(s) or other components. Note: In some lifts, a sudden quick movement at the beginning of the cycle is normal but should smooth out over the course of operation.
- Creating a side load on the actuator by pushing and/or pulling on the actuator, especially under load conditions, can cause damage to the actuator and or personal injury.
- Electronic components are designed to be plug-and-play type replacements. There is no need to access the internal mechanisms in these components. The actuator, battery, and all electronic components may only be disassembled by authorized factory-trained personnel. Disassembly of electronic components by improperly trained and unauthorized persons will void the warranty on those components.
- The operator of the pool lift should never allow the pool lift to hit any stationary object either vertically or horizontally while operating a pool lift, i.e. walls, ceiling, pool/spa bench, etc.
- For quick and accurate troubleshooting and inspections, we recommend having a voltage meter on hand as well as adequate test and replacement components for the unit.
- Remove the pool lift from service immediately if any unusual noise is heard coming from the actuator such as grinding, popping or sudden inconsistent movements.
- Immediately remove from service any lift that does not have 100% compliance with the inspection until the appropriate corrections have been made by a qualified person.

OVERALL VISUAL INSPECTION

Corrosion is common with anything in a swimming pool environment; even stainless steel. Proper and regular maintenance and cleaning will prevent most instances of corrosion. Replace any components that have corrosion caused pitting, cracking, thinning of the base metal, or other instances of corrosion that can't be eliminated and/or minimized through thorough cleaning. refer to product manual for additional maintenance instructions. Make comments and/or recommendations in the notes area below each section.

INSPECTION ITEM	PASS	FAIL	N/A
1. Anchor bolts, concrete around anchors, etc... are in good condition. No cracks in concrete, missing hardware, and anchors are installed correctly.			
2. The fixing anchor is installed correctly and is in good condition.			
3. Powder coated components are clean, corrosion free, and in good overall condition. Clean and repair any missing powder coat finish. Apply a coat of wax on the cleaned and repaired finish.			
4. Welded joints show no signs of cracks or fatigue. Replace any non-compliant components.			
5. Hardware is present, secure, and not corroded (bolts, nuts, bushings, etc...).			
6. All wiring routing clips are present and in good condition.			
7. All product warning labels are present and legible.			
8. Plastic seat, foot rest, and other plastic components are present and in good shape.			
9. Armrest hand grips are present, in good condition, and are well secured to the arms.			
10. Electronic and lift covers are present and in good condition (Lift covers are highly recommended for outdoor installations).			
11. Inspect rotation gears for any damage or missing teeth. Inspect rotation motor gear and alignment to ensure proper rotation of the lift.			
12. Inspect lift slings, seat belts, and other similar stitched and webbed or weaved materials for damage, fatigue, and/or deterioration. Replace at first sign of any fatigue or damage.			

PHYSICAL POOL LIFT TEST

A physical test of the lift function is required when doing an annual inspection. NOTE: The lift should always be run through a completed cycle without a load prior to each use. The annual test should be ran both empty and with a typical load weight to test the lifts functions. Make comments and/or recommendations in the notes area below each section.

INSPECTION ITEM	PASS	FAIL	N/A
1. With all connections connected, operate the lift through a complete cycle without weight to confirm the lift is functioning correctly. Note any strange noises and any unusual operation.			
2. Verify the lift was installed correctly; is the lift free to travel in and out of the pool without hitting anything like a wall, rail, bench, etc...			
3. Rotate lift at least 90 degrees over the water and back again if applicable.			

ELECTRONICS INSPECTION

Most electronic issues are due to pool electrical connections and/or a defective battery. Always remove the battery and place on the charger when the lift is not in use. NOTE: Even if your lift is equipped with a solar charger, it is a good idea to place the battery on the wall charger periodically or during extended periods of non-use or storage. Before testing the battery, charge overnight. A voltage meter is required for this section and a known GOOD handset controller is helpful. An Aqua Creek electronics test light is also highly recommended (order item F-01TL). Make comments and/or recommendations in the notes area below each section.

INSPECTION ITEM	PASS	FAIL	N/A
1. Battery contact points are free of damage and/or corrosion. Dielectric grease is present on all contact points to prevent corrosion. ALWAYS STORE BATTERY ON THE CHARGER WHEN LIFT IS NOT IN USE!			
2. Battery level indicator (if equipped) on control box is functioning properly when battery is installed on lift.			
3. Using the voltage meter check the voltage on the battery. A good battery should read between 24v and 28v DC. Note: Vito batteries will require a test cord. (F-VTC002)			
4. Using the test light, confirm continuity through the control box by plugging the light into the control box in place of the actuator and press the buttons on the controller and the control box itself (if applicable). If the light does not light up with the handset, replace handset with a known good one and try to see if it lights up.			
5. Once the battery is fully charged up to 24 to 28 volts, cycle the lift (up, down, rotate) counting the cycles. If the battery is in good condition, the lift should easily complete at least 15 cycles without a load.			

BASIC TROUBLESHOOTING

MAKE NOTE OF THE COMPLAINT, SUCH AS LIFT IS NOT OPERATING AT ALL, LIFT WON'T RETURN TO DECK, ETC...

MAKE NOTE OF ANY NOISES THE LIFT MAY BE MAKING, SUCH AS GRINDING, CLICKING, BUZZING, BEEPING. TRY TO IDENTIFY WHERE THE NOISES ARE COMING FROM, I.E. THE ACTUATOR, BATTERY CONTROL BOX, ETC...

IF A 'BEEPING' NOISE IS COMING FROM THE BATTERY, CHANCES ARE THE BATTERY IS DEAD OR NOT CHARGED ENOUGH TO OPERATE THE LIFT. IT IS DESIGNED TO HAVE AN AUDIBLE ALARM WHEN IT REACHES BELOW A CERTAIN CHARGE LEVEL. CHECK THE BATTERY VOLTAGE (SHOULD READ >24 VOLTS DC) AND EITHER CHARGE THE BATTERY OR REPLACE IT WITH A NEW ONE.

- A. LAY THE BATTERY ON A FLAT, DRY SURFACE.
- B. MAKE SURE YOUR VOLTMETER IS ON AND SET TO DC MODE.
- C. PUT ONE PROBE ON THE LEFT METAL PRONG AND ONE OF THE RIGHT METAL PRONG.
- D. THE VOLTMETER SHOULD READ 24 VOLTS OR HIGHER FOR A GOOD BATTERY. IF IT READS NEGATIVE IT JUST MEANS THE PROBES ARE ON BACKWARDS.
- E. SOMETIMES IT WILL TEST GOOD ON VOLTAGE BUT THE AMPERAGE IS INSUFFICIENT.
- F. THE BATTERY CAN BE TAKEN TO A BATTERY STORE TO HAVE AN 8 AMP LOAD TEST DONE TO CHECK FOR PROPER AMPERAGE.

NOTE: THIS CONTROL BOX COMES WITH AN "EMERGENCY STOP" BUTTON ON THE TOP. LIFT WILL NOT WORK WHILE THE BUTTON IS COMPRESSED. PLEASE ENSURE BUTTON IS NOT ENGAGED WHILE LIFT IS IN OPERATION. TO DISENGAGE, PRESS AND TWIST BUTTON UNTIL IT POPS UP.

1. IF THE BATTERY JUST CAME OFF THE CHARGER AND IS FULLY CHARGED, CHECK TO MAKE SURE THE CHARGER IS PLUGGED IN AND CHARGING PROPERLY. THIS CAN BE DONE WITH THE VOLTMETER.
 - G. VERIFY THE CHARGER IS PLUGGED INTO A WORKING WALL OUTLET.
 - H. PLACE ONE PROBE FROM YOUR VOLTMETER ON THE LEFT METAL PRONG ON THE CHARGER AND THE OTHER PROBE ON THE RIGHT METAL PRONG OF THE CHARGER. VOLTAGE SHOULD BE READING AT LEAST 27 VOLTS DC. IF IT READS NEGATIVE IT JUST MEANS THE PRONGS WERE PLACED BACKWARDS.
1. VERIFY THE CORDS ARE PROPERLY SEATED IN THE CONTROL BOX PORTS AND ARE NOT DAMAGED OR CORRODED. WE SUGGEST UNPLUGGING AND THEN PLUG BACK IN TO ENSURE THEY ARE SEATED PROPERLY. CLEAN ANY CORROSION OFF OF THE BATTERY TABS AND CORD PLUGS IF THERE IS ANY.

ONCE THESE STEPS ARE COMPLETE CALL AQUA CREEK PRODUCTS AND TALK TO TECHNICAL SUPPORT FOR FURTHER ASSISTANCE OR GO ONLINE AND WATCH OUR TROUBLESHOOTING VIDEOS: (888)687-3552
[HTTPS://AQUACREEKPRODUCTS.COM/HELP-CENTER/TROUBLESHOOTING/](https://aquacreekproducts.com/help-center/troubleshooting/)



LIMITED FIVE (5) YEAR WARRANTY

Aqua Creek Products, LLC (a.k.a. Aqua Creek) warrants to the original end user purchaser that products manufactured by Aqua Creek, when properly installed in accordance with assembly and installation instructions, and properly used and maintained, shall be free from defects in material and workmanship for a period of *five (5) years from the date of original purchase*, provided that Aqua Creek receives prompt notice in writing of any defect or failure and satisfactory proof thereof, with the following exception(s):

Exceptions:

All provided electrical components, including the linear actuator shall have the following warranty period:

Year 1-2: 100% Coverage

Year 3: 60% Coverage (Customer is responsible for 40% of replacement cost)

Year 4: 50% Coverage (Customer is responsible for 50% of replacement cost)

Year 5: 40% Coverage (Customer is responsible for 60% of replacement cost)

Hydraulic actuators and mesh slings shall have a warranty period of *one (1) year from the date of original purchase*.

Powder coat finish scratches, scrapes, corrosion, or dents from customers normal use, negligence, or abuse.

This warranty specifically excludes reimbursement for labor to remove, repair, or install the product and any shipping charges. These warranties do not cover any damages due to accident, force majeure, act of nature, misuse, abuse, negligence or failure to properly maintain any products, or normal wear and tear from day to day operations. In the event that any products are altered, repaired, or improperly installed or improperly used by anyone without the prior written approval by Aqua Creek, all warranties are void. **IMPORTANT: AMOUNT OF WEIGHT PLACED ON LIFT SHALL NOT EXCEED THE RATED LIFTING CAPACITY FOR THE LIFT. NEVER OPERATE THE LIFT UNDER LOAD IN A DRY POOL (WITH NO WATER IN THE POOL).** It is the responsibility of the lift owner to verify the weight of the patron for warranty claim purposes and to ensure that the lift is not overloaded. Non-payment for product to Aqua Creek may void warranty.

To initiate a warranty claim, the owner of an Aqua Creek product must provide the place of purchase, in writing, with a full description of the product, its serial number, the dates of purchase and installation, and the exact nature of the defect. Within thirty (30) days after receipt of a written warranty claim by Aqua Creek, and barring any unforeseen delays, the place of purchase will be notified of Aqua Creek's decision regarding the claim.

If requested by Aqua Creek, any defective product must be returned, freight prepaid by customer, to Aqua Creek's designated factory location or duly appointed distributor for inspection and/or repair. Aqua Creek will, at its option, repair or replace the failed or defective item, and deliver the repaired product or replacement to the customer. All freight charges are the responsibility of the customer. Products returned to Aqua Creek for which Aqua Creek provides replacement under this limited warranty shall become the property of Aqua Creek. A new warranty period shall NOT be established for the repaired or replaced products. Such products shall remain under warranty only for the remainder of the original warranty period on the original products purchased.

This written limited warranty constitutes the final, complete and exclusive statement of warranty terms. No person or organization is authorized to make any other specific or implied warranties or representations on behalf of Aqua Creek.

THE WARRANTIES SET FORTH HEREIN ARE IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE.

THE SOLE AND EXCLUSIVE REMEDIES FOR BREACH OF ANY AND ALL WARRANTIES WITH RESPECT TO THE PRODUCTS SHALL BE LIMITED TO REPAIR OR REPLACEMENT AT AQUA CREEK'S DESIGNATED FACTORY LOCATION, OR DULY APPOINTED DISTRIBUTOR, OR IN PLACE AT AQUA CREEK'S OPTION. IN NO EVENT SHALL AQUA CREEK'S LIABILITY EXCEED THE ENTIRE AMOUNT PAID TO AQUA CREEK BY THE ORIGINAL PURCHASER FOR THE FAILED OR DEFECTIVE PRODUCT. IN NO EVENT SHALL AQUA CREEK PRODUCTS, LLC BE LIABLE FOR ANY INCIDENTAL, CONSEQUENTIAL, SPECIAL, INDIRECT, PUNITIVE OR EXEMPLARY DAMAGES OR LOST PROFITS FROM ANY BREACH OF THIS LIMITED WARRANTY OR OTHERWISE.

THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS, WHICH MAY VARY FROM STATE TO STATE. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL, SPECIAL OR CONSEQUENTIAL DAMAGES, SO SOME OF THE ABOVE LIMITATIONS OR EXCLUSIONS MAY NOT APPLY TO YOU.